

PICARRO

Calibration strategies for cavity ring-down spectroscopy

Renato Winkler, Ph.D
Applications Scientist

HIGHGAS Stakeholder Workshop
November 14th, 2014

Greenhouse Gas instruments



CRDS Analyzer CO₂ CH₄ H₂O
Measurements in Air
Model G2301

The world's highest precision analyzer for the top three greenhouse gases

- Global #1 in precision and accuracy, and portability
- The lowest guaranteed drift of any instrument
- Unique water correction feature automatically reports dry mol fraction
- Innovative software featuring intuitive user interface & customization tools
- World class customer service and technical support




CRDS Analyzer
CO₂ + CO + CH₄ + H₂O
Model G2401

The world's only field-deployable analyzer capable of measuring the four main atmospheric trace gases simultaneously and continuously.


- Global #1 in precision, accuracy, and portability
- Capable of meeting WMO Data Quality Objectives for CO, CO₂ & CH₄
- Guaranteed lowest drift of any continuous greenhouse gas measurement instrument
- Unique water correction automatically reports dry gas mol fractions

GHG Measurements: The Mauna Loa CO₂ Record




The measurement cell – small and robust

← 30 cm →







Cavity Ring-Down Spectroscopy – Hardware

- Laser light at a wavelength that the molecule does absorb light





Stronger light absorbance -> shorter ring-down time



Pressure and temperature stability

Tightly control the conditions inside the cavity, so that variations of the environmental conditions* have no significant effect on:

- Cavity Temperature Stability ($\pm 0.007\text{ }^\circ\text{C}$)
- Cavity Pressure Stability ($\pm 0.00002\text{ atm}$)



*The analyzer is placed in an environmental chamber where the temperature is ramped up (then down), while measuring a cylinder of compressed air.



What Makes a Great GHG Instrument?

- **Stability, Stability, Stability!**
 - The more unstable the instrument is, the more you have to calibrate, and calibration gases are:
 - Costly
 - Difficult to deploy remotely (it's not so easy to ship gas bottles from NOAA to Greenland or Borneo or wherever the station is)
- **Small sample size**
 - The ability to measure a small amount of gas means that the expensive calibration standard lasts a long time

Cavity Ring-Down Spectroscopy – Low Drift

The long-term effects of this control and stability are clear

**Two Analytes,
One Month
Zero Calibration**




Picarro G2401 meets WMO /GAW specifications

Performance Specifications, in dry air	CO ₂	CO	CH ₄	H ₂ O
Precision (5 sec / 5 min, 1-σ) <i>Guaranteed for below specified range & conditions - reference gas not needed</i>	<150 ppb / <50 ppb	<30 ppb / < 2 ppb	< 1 ppb / < 0.7 ppb	< 200 ppm / < 50 ppm
Max Drift at STP (over 24 hrs / 1 month) (peak-to-peak, 50-minute average) <i>Guaranteed for below specified range & conditions - reference gas not needed</i>	150 ppb / 500 ppb	15 ppb / 50 ppb	1 ppb / 3 ppb	100 ppm \pm 5% of reading
Max Uncertainty Using Reference Gas* (1hr average, 2-sigma) <i>Meets WMO Data Quality Objective for GAW Stations</i>	< 50 ppb	< 2 ppb	< 1 ppb	n/a
Automated Determination of Dry Mol Fraction	Included	Included	Included	n/a
Operating Range	0 – 1000 ppm	0 – 5 ppm	0 – 20 ppm	0 - 7 %v H ₂ O / 39 °C dew pt (non-condensing)
Guaranteed Specifications Range	300 – 500 ppm	0 – 1 ppm	1 – 3 ppm	0 - 3 %v H ₂ O / 25 °C dew pt (non-condensing)
Measurement Interval	< 5 seconds	< 5 seconds	< 5 seconds	< 5 seconds
Rise/Fall time (10 - 90 % / 90 - 10%)	< 5 seconds	< 5 seconds	< 5 seconds	< 5 seconds

“Calibrate-ability” (that is, drift) in a CRDS Instrument

**(Actually, what follow is true of most optical
instruments)**

Calibration, is it needed?...yes it is




- 1)
Even a perfect instrument with a 1-sigma standard deviation of ~ 0
Does not report the “**Truth**” and has a non-zero **Bias**
- 2)
At Picarro we build instruments that have very high **precision** and very low **drift**

Precision:

Precision is the closeness of agreement between independent measurements of a quantity under the same conditions. It is a measure of how well a measurement can be made **without reference to a theoretical or true value.**


www.nde-org

How Calibration Is Performed




What is the physical basis of the difference?

The basic elements of optical spectroscopy



The basic elements of optical spectroscopy



What if the Spectrometer Isn't Perfect?

True concentration

Measured peak height (and proportionality k)

$$c = xkm + \varepsilon$$

Span ($x \sim 1$)



Offset ($\varepsilon \sim 0$)

- “SPAN” drift

- Incorrect temperature
- Incorrect Pressure
- Wavelength error
- Spectral interference from other species

- “OFFSET” drift

- “micro-etalons”, or baseline ripple
- Spectral interference from other species



Calibration of GHG Concentration Instruments

Drift in a G2401 (CFKADS)

True concentration


Measured peak height (and proportionality k)

$$c = xkm + \varepsilon$$

Span ($x \sim 1$)

Offset ($\varepsilon \sim 0$)

- “SPAN” drift dominates CH_4 and CO_2
 - Time scale is minutes to hours
- “OFFSET” drift dominates CO and H_2O (only at low water levels)
 - Time scale is hours to days



Typical Calibration Strategy for GHGs

- Calibrate with 3-5 bottles infrequently (start of life, yearly, etc) to verify instrument linearity over full range of measurements
- Calibrate in the field on a daily / weekly / monthly basis with 2-3 bottles
- “Target Tank” - Always have one known bottle that does not participate as a calibration reference to validate that the sampling system and calibration methodology gives the right answer

Calibration of $\delta^{13}\text{CH}_4$ and $\delta^{13}\text{CO}_2$ for atmospheric measurements


Testing the Instrument: Precision and Drift

Precision, $\delta^{13}\text{C}$ in CH_4
($1-\sigma$, 1 hr window)

< 0.8 ‰ guaranteed precision at > 1.8 ppm 5 min.
average

average

< 0.5 ‰ guaranteed precision at > 1.8 ppm, 15 minute
average



How Does Drift in an Infrared Isotope Spectrometer Manifest?

$$\delta_{13} = 1000 \left(\frac{r_{sample}}{r_{VPDB}} - 1 \right) \quad r_{sample} = c_{13}/c_{12} \quad c = xkm + \varepsilon$$

$$\delta_{true} = 1000 \left(\frac{(x_{13}k_{13}m_{13} + \varepsilon_{13})/(x_{12}k_{12}m_{12} + \varepsilon_{12})}{r_{VPDB}} - 1 \right)$$

Ignore ε_{12} when $m_{13} \ll m_{12}$, and $x_{13} \sim x_{12} \sim 1$

$$\delta_{true} = 1000 \left(\begin{array}{l} \frac{k_{13}m_{13}/k_{12}m_{12}}{r_{VPDB}} - 1 \\ + \left(\frac{x_{13}}{x_{12}} - 1 \right) \frac{k_{13}m_{13}/k_{12}m_{12}}{r_{VPDB}} \\ + \varepsilon_{13}/x_{12}k_{12}m_{12}r_{VPDB} \end{array} \right)$$

δ_{raw}

Delta span drift

[concentration]⁻¹ drift

$$\delta_{true} = \delta_{raw} + B(\delta_{raw} + 1000) + \frac{A}{^{12}\text{CH}_4}$$

A HI – LO Calibration Scheme

$$\delta_{true} = \delta_{raw} + B(\delta_{raw} + 1000) + \frac{A}{^{12}CH_4}$$


- *delta span drift* affects high and low concentration equally
- *conc⁻¹ drift* affects low concentration bottles more
- Calibration Strategy: Use two known bottles, a HI (C_{HI} and δ_{HI}) and a LO (C_{LO} and δ_{LO}), to determine A_i and B_i in each calibration period τ_i

Define: $\Delta_H \equiv \delta_{HI} - \delta_{HI\ raw}$ and $\Delta_L \equiv \delta_{LO} - \delta_{LO\ raw}$


$$A_i = \left[\frac{\Delta_H(\delta_{LO\ raw} + 1000) - \Delta_L(\delta_{HI\ raw} + 1000)}{\frac{(\delta_{LO\ raw} + 1000)}{C_H} - \frac{(\delta_{HI\ raw} + 1000)}{C_L}} \right]_i$$

$$B_i = [(\Delta_H - A_i/C_H)/(\delta_{HI\ raw} + 1000)]_i$$

Testing the Scheme – 40 days and 40 nights




- Hourly measurement cycle
 - HI = 10.07 ppm / -41.7 ‰ (10 min)
 - LO = 1.80 ppm / -44.5 ‰ (25 min)
 - TARGET = 2.03 ppm (25 min)




Calibration – [conc]⁻¹ Drift Dominates

$$\delta_{true} = \delta_{raw} + B(\delta_{raw} + 1000) + \frac{A}{CH_4}$$




>90% of the drift can be corrected using a single bottle and a $A / ^{12}CH_4$ correction factor

Target Tank Stability After Hourly Calibration



$\sigma_{1 \text{ hr}} = 0.204 \%$
 $\sigma_{24 \text{ hr}} = 0.046 \%$



Recommended Calibration Strategy for isotopic measurements

- Calibrate with 3 (or more) ref. gases infrequently (start of life, yearly, etc) to verify instrument linearity over full range of isotopic and concentration range, and to characterize your working standards
- Calibrate in the field on every couple of hours / 2-times daily with 2 working standards bracketing the range in delta values and concentration (HI LO scheme)
- “Daily working standards” - Always have 1-2 bottles that has an isotopic composition close to the samples, and acts as a quality control (daily offset and drift)

Thank you for your attention

Questions?

- Feel free to contact me:

rwinkler@picarro.com

Picarro analyzers: used and recommended by the world's leading scientists, in industry and academia, and by governmental bodies*.



*selected Picarro customers